




# BND Inc.

Intertek  
ISO:13485:2016 - MDSAP

Stretch Your Healthcare Dollar

# COVID-19 and Flu Viral Collection Kits



## Virus Sampling Kit

By Yacon Biology

- 3.5 ml of Hanks balanced salt solution based viral transport media
- Self-standing 10ml Tube with Screw Cap
- Sterilized Kit contains Tube with Medium & Swab
- Head material: Nylon flocked swabs
- Handle material: ABS
- Breakpoint distance: 80mm or 100mm
- Total length: 150 mm
- Storage: Store at room temperature (5-25°C) away from light.
- Shelf Life: 12 months
- Nasal or oral swab kits
- 20 Kits/Box

Description	Item
Nasal Virus Sampling Kit	STC-110
Oral Virus Sampling Kit	STC-120

## Viral Transport Medium Kit

By Miraclean Technology Co Ltd.

- 3.0 ml of Hanks balanced salt solution based viral transport media
- Self-standing 10ml Tube with Screw Cap
- Sterilized Kit contains Tube with Medium & Swab
- Head material: Nylon flocked swabs
- Handle material: ABS
- Breakpoint distance: 82 or 95 mm
- Total length: 150 mm
- Storage: Store at room temperature (5-25°C) away from light.
- Shelf Life: 24 months
- Nasal or oral
- 50 kits/pack for transport medium kits or 100 swabs/pack for standalone swabs



Description	Item
Nasal Viral Transport Medium Kit	STC-100
Standalone flocked Swab	STC-200
Oral Viral Transport Medium Kit	STC-300
Children's Standalone Swabs	STC-400

## Saliva Collection Kits and Tubes

- Collection tubes in multiple formats
- Collection tubes with or without transport media
- Collection funnels in varying sizes and configurations
- Sterile and individually packaged




## Biohazard Specimen Transport Bag

Technologists Choice™ Biohazard Specimen Transport Bag

- Re-closable, avoids contamination of inner packed items
- Colour, size, print and packaging can be customized
- Optional document compartment with flap
- Single zipper, 2 mil thickness
- 500 Pieces/Case

Description	Item
Clear with red and black print 6 x 9 x 2 mil	DIS-028



**ORDER TODAY!**

1-800-668-4033 • CustomerService@bndinc.com